Geratherm[®]
smart

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CE 0197

**Digital Blood Pressure Monitor
for the Upper Arm**

Geratherm[®]
smart

EN

INSTRUCTION MANUAL

GT-1775

CE0197

Contents

General description	41
Warranty	42
Measurement principle	42
Precautions	47
Important information for the user	48
Description of the device	50
Information about the device	50
Explanation of the display	51
Preparation for use	52
Inserting and replacing the batteries	52
Setting the time and date	53
Attaching the cuff	56
Posture during measurement	57
Functions	58
Measuring blood pressure	58
Detecting an irregular heartbeat	59
Retrieving stored readings	61
Deleting stored readings	62
Care and maintenance	63
Error messages	64
Technical data	66
Quality standards	68
Index of symbols	68
Information on electromagnetic compatibility (EMC)	69
Appendix	70

EN

General description

Thank you for choosing the Geratherm® smart blood pressure monitor.

Before using your device for the first time, please read this instruction manual carefully.

This blood pressure monitor uses the oscillometric method of measuring the systolic and diastolic blood pressure and the heart rate. The measurements are taken at the upper arm and all values appear on an LCD display screen.

This blood pressure monitor complies with the standards for medical devices. The blood pressure monitor can be used by adults with an arm circumference of 22 cm to 42 cm.

 Before using this device, pregnant women must consult a physician and the suitability of this measuring method must be assessed. In some cases, the measurements may deviate. Pregnant women are therefore recommended to have control measurements conducted by their physician, since clinical validation has not been conducted with pregnant women.

Persons with implanted electrical devices such as heart pacemakers and defibrillators must consult a physician regarding the suitability of this measuring method.

EN

Contents:

1. Blood pressure monitor Geratherm® smart (GT-1775)
2. Cuff (22 cm - 32 cm)
3. 4 x type AAA batteries (1.5 V DC)
4. Instruction manual
5. storage bag

Note:

A 22 - 42 cm size cuff for the Geratherm® smart can also be purchased from the dealer.

Warranty

This blood pressure monitor is backed by a warranty against any manufacturing defects for three years from the date of purchase, subject to normal usage. If your blood pressure monitor does not function properly due to defective parts or assembly, we will repair it free of charge. Except for the cuff, all parts of the blood pressure monitor are covered by this warranty. Damage caused by improper use or handling of this blood pressure monitor is not covered by the warranty.

Measurement principle

EN

This product uses the oscillometric measuring method to determine blood pressure. Before every measurement, the unit establishes a “zero pressure” equivalent to the air pressure. Then it starts inflating the cuff. As it does so, the unit detects pressure oscillations generated by the pulsation of the blood flow, thus determining the systolic and diastolic pressure, as well as the heart rate.

Precautions

- The device is intended for private use only.
- These instructions and the device are no substitute for visiting your doctor. Neither the information contained here nor this device may be used to diagnose or treat health problems or to prescribe medicines. If you have a medical problem or suspect that you may have such a problem, please seek immediate advice from your doctor.
- This device is not intended for infants, pregnant women, patients with implanted electronic devices, patients with pre-eclampsia, premature ventricular beats, atrial fibrillation, peripheral artery disease and patients who are undergoing intravascular therapy or who have an arteriovenous shunt, or persons who have had a mastectomy. A doctor must be consulted prior to using the device on any of the persons listed above, and the suitability of the measuring method must be assessed. Geratherm Medical AG will not accept liability for any loss or damage that may occur.
- The device is not intended for measuring blood pressure in young children. Consult your doctor prior to using the device on older children (minimum arm circumference 22cm).
- The device is not intended for the transportation of patients outside medical establishments.
- This device is intended for the non-invasive measuring and monitoring of arterial blood pressure. It is not intended for use on other extremities than the arm or for any other functions other than measuring blood pressure.
- This device is not suitable for continuous monitoring during medical emergencies and surgery, since the circulation of blood in the patient's arms and fingers will be disrupted owing to a lack of blood. They may become numb, swollen and even blueish-red.

EN

Precautions

- Do not mistake self-monitoring for self-diagnosis. This device enables you to monitor your blood pressure. Do not undertake any therapeutic measures on the basis of measurements you have taken yourself.
- Under no circumstances should you change the dosage of the medication that your doctor has prescribed for you. Consult your doctor if you have any questions about your blood pressure.
- If you take any medication, ask your doctor about the best time for measuring your blood pressure.
- In patients with frequent heart rate irregularities, such as premature atrial or ventricular contractions or atrial fibrillation, the measuring results may be inaccurate. Please consult your doctor concerning the results.
- Avoid touching the patient while inserting or changing the batteries.
- Wait 30 to 45 minutes before measuring if you have recently consumed a drink containing caffeine or have smoked a cigarette.
- Relax for at least 5 to 10 minutes before each measurement.
- Keep the device out of the reach of infants, children and pets. The inhalation or swallowing of small parts is dangerous and can be fatal.
- In order to avoid accidental strangulation by the tubing of the cuff, keep the product away from children.
- During measurement, please ensure that the blood circulation is not hindered.
- If you experience any discomfort during measurement, such as pain in the arm or other complaints, press the "START/STOP" button to deflate the cuff. Release the cuff and remove it from your arm.

Precautions

- In the rare event of a fault causing the cuff to remain fully inflated (cuff pressure >300 mmHg) during measurement, open the cuff immediately. Prolonged high pressure applied to the arm may result in an ecchymosis.
- Too frequent and consecutive measurements can lead to disturbances in the blood circulation and cause injury.
- Do not wrap the cuff around an arm in which an arteriovenous shunt has been implanted or to which other medical devices are attached, as this could lead to temporary functional failure of the medical devices being used.
- Avoid any kinking, crushing or constricting of the connecting tube of the cuff during use, as this could lead to a continuous increase in the cuff pressure, resulting in interruption of the blood flow and causing injury to the patient.
- In order to avoid inaccurate measurements, do not expose the blood pressure monitor to strong electromagnetic fields, stray radiation signals or rapid transient electrical disturbances/burst transmissions.
- **Important!**
Do not use the cuff on damaged skin.
- Wireless communications devices such as WiFi, mobile phones, cordless phones and their base units or walkie-talkies, may impair the functioning of the blood pressure monitor. Therefore, a minimum separation distance should be maintained. Please refer to Table 4/9 in the appendix.
- Please do not use this device if you are allergic to polyester, nylon or synthetic materials.
- Do not wash the cuff in a machine! Wash only by hand!

EN

Precautions

- The service life of the cuff can vary according to the frequency of washing, the condition of the skin and the manner of storage. The typical service life is approximately 10,000 measurements.
- The materials used in the cuff have been tested and have been found to comply with the requirements of the standards ISO 10993-5:2009 and ISO 10993-10:2010. They do not have the potential to cause any contact sensitisation or irritation.
- Do not connect the air tube to any other medical equipment, as this can lead to high pressure or conduct air into intravascular systems, which may lead to injury.
- This blood pressure monitor is not designed for use in a type AP or APG gas-filled atmosphere or in the presence of inflammable anaesthetic gases with air or oxygen or nitrogen oxide.

- Before each use, check that the device is undamaged and is in the correct operating mode. Do not use the unit if it is damaged! Use of a damaged device can cause injury or give inaccurate readings which may constitute a serious hazard.

EN

- Allow the blood pressure monitor to warm up or cool down for at least 30 minutes if it has been stored at minimum or maximum storage temperature. Only when operating temperature has been reached is the unit ready for use.
- Please use the device under the environmental conditions described in the instructions. Otherwise the performance and service life of the blood pressure monitor will be adversely affected.
- When it is not in use, please store the blood pressure monitor in a dry room and protect it from extreme humidity, heat, fluff, dust and direct sunlight. Do not place any heavy objects on top of the product packaging. Also remove the batteries.

Precautions

- The blood pressure monitor contains sensitive components and must be treated with care. Therefore, always comply with the storage and operating conditions described above.

- In order to test the measuring accuracy of the Geratherm® smart when necessary, please contact the manufacturer or your dealer.
- Please use only accessories that have been approved by the manufacturer. Otherwise, damage may be caused to the device, the user may suffer injury, or inaccurate measurements may occur.
- Dispose of the unit, accessories and components in accordance with local regulations.
- In the event of malfunction, do not attempt to repair the device yourself. Only permit repairs to be carried out by authorised technicians.
- Never attempt to open or repair the unit while it is in operation.
- In the event of an unexpected error message or problems with the product, please contact your dealer or Geratherm Medical AG.
- Please notify Geratherm Medical AG of any unexpected occurrences or incidents.
- The Geratherm® smart blood pressure monitor must be used only for its intended purpose of measuring blood pressure. The manufacturer cannot be held liable for any damage caused by improper use or failure to comply with the safety notes given in the instruction manual.

EN

Important information for the user

What is blood pressure?

One force is created by the heart as the chamber of the heart, called the ventricle, forcibly ejects blood into the blood vessels and through the vascular system. Another force is created by the arteries as they resist the blood flow. Blood pressure is the result of these two forces.

Is my blood pressure normal?

To evaluate your blood pressure, see the following blood pressure classification chart published by the WHO (World Health Organization).

Please note:

Deviations are very rare. Please talk to your doctor if the results you measure are consistently outside the normal range.

Blood pressure classification	Systolic mmHg	Diastolic mmHg
Optimal	< 120	< 80
Normal	120 - 129	80 - 84
High - Normal	130 - 139	85 - 89
Grade 2 Hypertension	140 - 159	90 - 99
Grade 1 Hypertension	160 - 179	100 - 109
Grade 3 Hypertension	≥ 180	≥ 110

EN

What is meant by systolic and diastolic pressure?

Systolic pressure is the higher value, measured at the height of the heart's contraction. Diastolic pressure is the lower value, representing the pressure when the heart relaxes.

What about low blood pressure?


Generally speaking, a lower blood pressure reading is better as long as it does not cause symptoms such as dizziness or fainting.

Important information for the user

Why does my blood pressure fluctuate over the course of the day?

An individual's blood pressure varies throughout the day. It is influenced by a wide range of factors. Therefore, always measure your blood pressure under the same conditions.

The following chart shows possible blood pressure fluctuations over a 24-hour period.

The following factors will influence blood pressure readings and cause variations:

- Bathing
- Drinking alcohol
- Moving
- Eating
- Thoughts
- Conversation
- Exercise
- Stress
- Temperature changes
- Smoking etc.


Does it matter whether the measurement is taken on the left or the right arm?

Measurements can be taken on the left or right arm. However, further measurements should always be taken on the same arm - preferably on the arm on which higher blood pressure is measured.

EN


Description of the device

Information about the device

List of components of the blood pressure measuring system

- 1 Cuff
- 2 Air tube
- 3 Electronic module
- 4 Pump
- 5 Valve


Note:


Only use the cuffs supplied by Geratherm® specifically for this device. A cuff designed for an arm circumference of 22 - 42 cm is also available as an option.

Description of the device

Explanation of the display

Symbol	Description	Explanation
SYS	Systolic pressure	Higher BP value in mmHg
DIA	Diastolic pressure	Lower BP value in mmHg
PUL/min	Pulse rate	Pulse in beats per minute
mmHg	mmHg	Unit of measurement for BP (1 mmHg = 0.133 kPa)
	Blood pressure indicator (level)	Indicates the BP level
	Weak battery	Batteries are weak and must be replaced.
	Irregular heartbeat	The device has detected an irregular heart beat during measurement.
	Current time	Day/Month/Year Hour : minute
	Pulse	The blood pressure monitor records the pulse rate (pulses per minute) during measurement.
AVG	Average value	Average Blood pressure readings

EN

Preparation for use

Inserting and replacing the batteries

Battery operation


4 x type AAA batteries (1.5 V DC)

Please note:


For the greatest efficiency, and in order to protect your device, please use only the batteries specified.

Inserting and replacing the batteries

- Open the battery compartment cover.
- Insert the batteries as shown, ensuring correct polarity.
- Close the battery compartment.

Change the batteries in the event of the following:

-  is displayed.
- The display screen is dim.
- The display screen fails to light up.

EN

Please note:


- Do not mix new and old batteries. Do not mix different types of batteries.
- Do not dispose of batteries in fire. The batteries may explode or leak.
- Remove the batteries if the device is expected to be out of service for some time.
- Empty batteries are damaging to the environment. Do not dispose of them in household waste.
- Dispose of batteries in accordance with your local recycling regulations.

Preparation for use

Setting the time and date


It is important to set the clock before using your blood pressure monitor, so that a time stamp can be assigned to each reading that is stored in the memory. (Setting range for the years 2017 - 2057; Time format: 24 hours).

1. With the device switched off, hold down the “START/STOP” button for 3 seconds in order to enter the mode for setting the year.


EN

2. Press the “M” button to change the [YEAR]. Each time you press, the number will increase by one.


Preparation for use

3. Press the "START/STOP" button in order to confirm the [YEAR]. The device will then switch to the mode for setting the [MONTH] and [DAY].


4. Repeat steps 2 and 3 in order to set the [MONTH] and [DAY].

EN


Preparation for use

- Repeat steps 2 and 3 in order to set the [HOUR] and [MINUTE].

- After the hours and minutes have been set, the LCD screen will then first display the word "DONE" and then the device will switch off.


EN

Preparation for use

Attaching the cuff

1. Before taking a blood pressure measurement, remove any jewellery, such as watches and bracelets, on the arm that is to be used for measuring.

Note:

If your doctor has diagnosed poor circulation in your left arm, use the right arm.


2. In order to be able to wrap the cuff around the arm, push or roll your sleeve up to expose the bare skin. Make sure your sleeve is not too tight. Hold your arm with the palm of your hand facing upwards and attached the cuff to the upper arm over the main artery. The air tube must be on the inner side of the arm in line with the little finger. On the cuff there is also a printed label that you can use for guidance.

EN

Note:

You can find the main artery by pressing two fingers approximately 2 centimetres above the bend of the elbow on the inner side of the arm. The point where the pulse can be strongly felt is the location of the main artery.

3. The cuff should be fitted snugly comfortably but should not constrict the arm. It should be possible to insert a finger between the cuff and the arm.

Preparation for use

Posture during measurement

1. Sit comfortably and lay your arm with the cuff on a flat surface.

Lay your elbow on the table so that the cuff is level with your heart.

The palm of your hand should face upwards. Sit upright and take 5 or 6 deep breaths.

Helpful tips for patients, particularly patients with hypertension:

- Rest for five minutes before measuring.
- Wait at least three minutes between measurements. This will enable your circulation to recover.
- Take the measurement in a quiet room.
- Try to relax during measurement; do not move or speak.
- Attach the cuff at a position level with your heart.
- Sit comfortably. Do not cross your legs but place your feet flat on the floor, and lean against the back of the chair.
- In order to ensure comparability of measurements, try to measure under similar conditions each time. For example, measure at the same time each day, using the same arm, or following the instructions of your doctor.


EN

Functions


Measuring blood pressure

1. With the device switched off, press the “START/ STOP” button to switch the device on. It will begin measuring.

LCD display


The value zero appears


EN

Pumping and measuring begin

Results are displayed and stored in the memory

2. Press the “START/STOP” button to switch the device off. **If you do not do this, the device will switch off automatically after one minute.**

Detecting an irregular heartbeat

An irregular heartbeat (arrhythmia) is detected if the cardiac rhythm fluctuates while the device is measuring systolic and diastolic blood pressure. During each measurement the device records the intervals between heartbeats and calculates the mean interval value. If there are major irregularities, the arrhythmia icon appears on the display screen.

'Major irregularities' are when two or more pulse intervals occur in which the difference between these intervals and the mean interval value deviates by more than $\pm 25\%$ from the mean interval value, or when there are four or more pulse intervals in which the difference between each interval and the mean interval value deviates by more than $\pm 15\%$ from the mean interval value. In these cases, the icon for arrhythmia will appear on the display screen.

PLEASE NOTE

The display icon “” indicates that a pulse irregularity consistent with an irregular heartbeat was detected during measurement. Usually this is NOT a cause for concern. However, if the symbol appears often, we recommend you seek medical advice. Please note that the device does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage.

Blood pressure indicator



The blood pressure indicator on the display screen shows which category your personal blood pressure falls into. A “P” on the display screen indicates optimal blood pressure.

Blood pressure (mmhg)	Level						
	optimal	normal	good	acceptable	moderate	severe	
SYS	< 120	120 - 129	130 - 139	140 - 159	160 - 179	≥ 180	
DIA	< 80	80 - 84	85 - 89	90 - 99	100 - 109	≥ 110	
							

Functions

Retrieving stored readings

1. With the device switched off, press the "M" button to view the average (AVG) of the last three stored values. If there are fewer than 3 groups of stored values, the latest stored value will be shown.
2. Press the "M" button again to retrieve a specific stored value.

The stored value will appear alternately with the date and time the reading was taken.

1 / 5

The sequential
no. is no. 1
5 stored values
in total

D M Y
1 / 1 / 17

The relevant
date is
1 January 2017

6:00

The relevant
time is
6 a.m.


EN

Functions

Deleting readings

You can delete all readings for the selected user by means of the following steps:

1. Press and hold the “M” button for 3 seconds when the device is already in memory mode. “DEL ALL” will then flash on the display screen.

2. Press and hold the “START/STOP” button to confirm the deletion process. The device will then display “dEL donE” and switch off automatically.


EN

Note:

To exit the Delete mode without deleting readings, press the “START/STOP” button again.

3. If you press the “M” button and there are not yet any readings in the memory, the display screen will appear as shown here:

Functions

Inaccurate readings may be obtained if the measurement is taken:

- within 1 hour after eating or drinking
- immediately after smoking
- within 20 minutes after having a bath
- when speaking or moving fingers
- in a cold environment
- when you have an urge to urinate

Care and maintenance


To keep your device in good working order, please comply with the following:

EN

- Store it in a dry place
- Protect it from direct sunlight, humidity and heat
- Avoid shaking or dropping the unit
- Use a dry cloth to remove dust and dirt
- Do not use any abrasive cleaners or solvents
- After contact with water, wipe immediately with a dry cloth
- Avoid dusty environments and fluctuating temperatures

Error messages

This section contains a list with frequently asked questions about problems and error messages that may occur with the blood pressure monitor. Please first check whether your problem can be solved with the aid of the tips given here. If the device still does not function as expected, please contact our customer service department.

Symptom	Display screen	Cause	Solution
No power	Display screen does not light up.	Batteries are empty.	Replace with new batteries.
		Batteries have been inserted incorrectly.	Insert the batteries correctly. Check the polarity! See section on "Inserting and replacing batteries"
Weak batteries	Display is dim or shows 	Batteries are weak.	Replace with new batteries.
Error message	E 01 appears	The cuff is too tight or too loose.	Re-fasten the cuff and measure again.
	E 02 appears	The device has detected movement during the measuring process.	Movement can influence the measurement. Relax for a moment and then measure again.

Error messages

Error message	E 03 appears	No pulse signal has been detected during measurement.	Loosen clothing around the arm and measure again.
	E 04 appears	Processing of the measurement has failed.	Relax for a moment and then measure again.
	EExx appears on the display	A calibration error has occurred. (XX May be a digital symbol such as 01, 02 etc. if the same situation results in the same calibration error)	Repeat the measurement. If the problem persists, contact your dealer or our customer service department for further support.
Warning message	"out" appears	Outside the measuring range	Relax for a moment. Re-fasten the cuff and measure again. If the problem persists, consult your doctor.

EN

Technical Data

Power supply	Battery operation: 4 x AAA batteries (1.5 V DC)
Display mode	Backlit LCD display screen approx. 65 mm x 50 mm
Measuring mode	Oscillographic testing mode
Measuring range	Rated cuff pressure: 0 mmHg ~ 299 mmHg Measurement pressure: SYS: 60 mmHg ~ 230 mmHg DIA: 40 mmHg ~ 130 mmHg Pulse rate: (40-199) beats per minute
Accuracy	Temperature: 5 °C - 40 °C to within ± 3 mmHg (0.4 kPa) Pulse value: ± 5 %
Operating conditions	A temperature range between: +5 °C and +40 °C A relative air humidity range of: 15 % to 90 %, non-condensing but not requiring water vapour partial pressure of more than 50 hPa An atmospheric pressure range of: 700 hPa to 1060 hPa
Storage and transportation conditions	Temperature: -20 °C to +60 °C A relative air humidity range of ≤ 93 %, non-condensing, with water vapour pressure up to 50 hPa
Cuff sizes	approx. 22 cm ~ 32 cm or 22 cm ~ 42 cm
Weight	approx. 225 g (without batteries and cuff)
Dimensions	approx. 120.2 mm x 108.2 mm x 68.5 mm

Technical Data

Mode of operation	Continuous operation
Degree of protection	Type BF applied part
Protection against ingress of water	IP21, this means: The device is protected against penetration by solid objects ≥ 12.5 mm in diameter and against the ingress of water in the form of vertically falling drops.
Device classification	Battery mode: medical electrical device with internal battery operation
Software version	A04

Important note:

No modification of this equipment is permitted.

EN

Quality standards

Geratherm Medical AG is certified in accordance with DIN EN ISO 13485 and fulfils the requirements of Council Directive 93/42/EEC. It is therefore authorised to apply the label CE0197 Notified Body: TÜV Rheinland LGA Products GmbH).

Index of symbols

EN

	Follow instruction manual		Type BF equipment
	Keep dry		Manufacturer
	Store at air humidity levels between 0 % and 93 % R.F.		(YYMMXXX; Year/Month/Serial number)
	Store in temperatures between -20 °C and +60 °C		Caution! Comply with safety regulations in the instruction manual
	The device must not be disposed of with household waste		Direct current

Information on electromagnetic compatibility (EMC)

Electronic devices such as PCs and mobile phones can lead to medical devices in operation being exposed to electromagnetic interference from other devices. This can lead to malfunction of the medical device and create a potentially unsafe situation.

Medical devices should also not interfere with any other equipment.

The EN 60601-1-2 standard regulates the requirements for EMC (electromagnetic compatibility) and defines the levels of immunity to electromagnetic interference and the maximum electromagnetic emission levels for medical devices.

This blood pressure monitor, which is manufactured by Geratherm Medical AG, complies with the EN 60601-1-2 standard in relation to both immunity and emissions.

Nevertheless, special precautions should be observed:

Please only use the device indoors and not in the vicinity of mobile phones or microwave ovens.

This device should not be used adjacent to or stacked with other devices; and if it has to be used adjacent to or stacked with other devices, it should be observed to verify normal operation in the configuration in which it is being used.

EN

Warning:

The use of accessories, transducers and cables other than those specified or provided by the manufacturer of this device may result in increased electromagnetic emissions or reduced electromagnetic immunity of the device and may cause it to malfunction.

Appendix

Guidance and manufacturer's declaration – electromagnetic emissions		
The device is intended for use in the electromagnetic environment specified below. The customer and/or user of the device should assure that it is used in such an environment.		
Emission test	Compatibility	Notes concerning the electromagnetic environment
RF emissions CISPR 11	Group 1	RF energy is used only to maintain the device's operation. Therefore, its RF emissions are so low that it is not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Compatible	

Guidance and manufacturer's declaration – electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer and/or user of the device should assure that it is used in such an environment.			
Immunity test	IEC 60601 - Test level	Compatibility level	Notes concerning the electromagnetic environment
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient / Burst IEC 61000-4-4	Power supply lines: ± 2 kV Input/output lines: ± 1 kV	Power supply lines: ± 2 kV	Mains power quality should be that of a typical commercial or hospital environment.
Surge voltage IEC 61000-4-5	Line(s) to line(s): ± 1 kV Line(s) to earth/ground: ± 2 kV 100 kHz repetition frequency	Line(s) to line(s): ± 1 kV 100 kHz repetition frequency	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % 1 cycle and 70 % 25/30 cycles Single-phase: at 0 0 % 300 cycles	0 % 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % 1 cycle and 70 % 25/30 cycles Single-phase: at 0 0 % 300 cycles	Mains power quality should be that of a typical commercial or hospital environment.
Power frequency (50 Hz / 60Hz) Magnetic field IEC 61000-4-8	30 A/m 50 Hz / 60 Hz	30 A/m 50 Hz / 60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Appendix

Recommended separation distances between portable and mobile RF communications equipment and the device or system.			
The device or system is intended for use in the specified electromagnetic environment, in which radiated RF disturbances are under control. The customer and/or user of the device can help prevent electromagnetic interference by keeping the device at a minimum distance from portable and mobile RF communications equipment (transmitters) as recommended below, according to the maximum output power of the communications equipment:			
Rated maximum output power of the transmitter in watts (W)	Separation distance according to the frequency of the transmitter (m)		
	150 kHz to 80 MHz $d = 3,5 \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \sqrt{P}$	800 MHz to 2,7 GHz $d = 2,3 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,37	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
For transmitters rated at a maximum power output not listed above, the recommended separation distance can be estimated using the equation in the relevant column, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1:	At 80 MHz and 800 MHz the higher frequency range applies.		
NOTE 2:	These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.		

Appendix

Guidance and manufacturer's declaration – electromagnetic immunity				
The device is intended for use in the electromagnetic environment specified below. The customer and/or user of the device should assure that it is used in such an environment.				
Immunity test	IEC 60601 – test level	Compatibility level	Notes concerning the electromagnetic environment	
Conducted RF IEC 61000-4-6	150 kHz to 80 MHz: 3 Vrms 6 Vrms (in ISM and amateur radio bands) 80 % Am at 1 kHz	150 kHz to 80 MHz: 3 Vrms (in ISM and amateur radio bands) 80 % Am at 1 kHz	Portable and RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation appropriate for the frequency of the transmitter. Recommended separation distance: $d = 0,5 \sqrt{P}$ $d = 1,2 \sqrt{P}$	
Radiated RF IEC 61000-4-3	10 V/m, 80 % Am at 1 kHz	10 V/m, 80 % Am at 1 kHz	80 MHz to 800 MHz: $d = 1,2 \sqrt{P}$ 800 MHz to 2,7 GHz: $d = 2,3 \sqrt{P}$	where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: <div style="text-align: right;">  </div>
NOTE 1:		At 80 MHz and 800 MHz the higher frequency range applies.		
NOTE 2:		These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.		

- a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.
- b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

EN

Guidance and manufacturer's declaration – electromagnetic interference immunity							
The device is intended for use in the electromagnetic environment specified below. The customer and/or user of the device should assure that it is used in such an environment.							
	Test frequency (MHz)	Band a) (MHz)	Service a)	Modulation b)	Modulation b) (W)	Distance (m)	Immunity test level (V/m)
Radiated RF IEC 61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to wireless RF communications equipment)	385	380-390	TETRA 400	Pulse modulation b) 18 Hz	1,8	0,3	27
	450	380-390	GMRS 460 FRS 460	FM c) \pm 5 kHz deviation 1 kHz sine 0.2	2	0,3	28
	710	704-787	LTE Band 13, 17	Pulse 2 Modulation b) 2 0.2 217 Hz	0,2	0,3	9
	745						
	780						
	810	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation b) 18 Hz	2	0,3	28
	870						
	930						
	1720	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT LTE Band 1, 3, 4, 25 UMTS	Pulse modulation b) 217 Hz	2	0,3	28
	1845						
	1970						
	2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation b) 217 Hz	2	0,3	28
	5240	5100-5800	WLAN 802.11 a/n	Pulse modulation b) 217 Hz	0,2	0,3	9
	5240						
	5785						

NOTE: If it is necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME device or ME system may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

- a) For some services, only the uplink frequencies are included.
- b) The carrier shall be modulated using a 50 % duty cycle square wave signal.
- c) As an alternative to FM modulation, 50% pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

The MANUFACTURER should consider reducing the minimum separation distance, based on RISK MANAGEMENT, and using higher IMMUNITY TEST LEVELS that are appropriate for the reduced minimum separation distance. Minimum separation distances for higher IMMUNITY TEST LEVELS shall be calculated using the following equation:

$$E = \frac{6}{d} \sqrt{P}$$

Where P is the maximum power in W, d is the minimum separation distance in m, and E is the IMMUNITY TEST LEVEL in V/m.

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